

# Efficient assays with orientation-specific binding of antibodies

Daphne Areskoug, Per Kjellin, Åsa Frostell, and Lena Vinterbäck

GE Healthcare Bio-Sciences AB, SE-751 84 Uppsala, Sweden

## Introduction

Many researchers are facing increasing challenges in assay development working with complex sample matrices, difficult targets, and often with the need to produce results at a faster pace. Antibody applications are supported by a number of Biacore™ sensor chips and reagents enabling rapid assay development. This poster presents three newly developed pre-immobilized sensor chips: Sensor Chip Protein A, Sensor Chip Protein G, and Sensor Chip Protein L. These sensor chips enable orientation-specific binding of antibodies for subsequent antibody concentration and other analyses and have a high binding capacity. The sensor chips use the same recombinant surface ligands as for GE's chromatography-based purification products.

Future Biacore consumables development follows a product category called Biacore Extend, with the purpose of development of sensor chips and reagents at a faster pace than previously for researchers working in nonregulated environments.



## Biacore Sensor Chip Protein A

The surface of Sensor Chip Protein A consists of a carboxymethylated dextran matrix with a recombinant Protein A variant covalently attached. Sensor Chip Protein A has the same surface ligand as is used for MabSelect SuRe™ protein purification products. This recombinant Protein A variant binds only to the heavy chain within the Fc region, which ensures a specific orientation of the antibody on the



surface. It may also be useful in antibody purification optimization since it does not bind to Fab fragments, as native Protein A does. Antibodies from several mammalian species bind to the surface, most notably human antibodies of the subclasses IgG<sub>1</sub>, IgG<sub>2</sub>, and IgG<sub>4</sub>. Sensor Chip Protein A is suitable for concentration analysis and characterization work involving kinetic analysis.

## Concentration analysis

### High sensitivity and wide dynamic range

Limit of detection (LOD), limit of quantitation (LOQ), and dynamic range were determined using a calibration curve for Avastin™, an IgG<sub>1</sub>-based mAb for the treatment of cancer, which inhibits the growth factor VEGF-A. A dynamic range from ng/ml to mg/ml was obtained (Fig 1). The colored lines show concentrations used for creation of the concentration curves. Colored areas of the same color as the lines indicate good fit to the calibration curve (< 10% deviation). LOD and LOQ were determined at 30 and 180 s injection times; both were improved at 180 s injection time (Fig 1).

Injection time (s)	30	180
LOD (ng/mL)	1.0	0.25
LOQ (ng/mL)	8.0	2.0




Fig 1. Avastin captured on Sensor Chip Protein A, response units as a function of Avastin concentration.

### High reproducibility

A sample series from two fed-batch cultures of a mAb-producing CHO cell line were analyzed. Analyses were performed on two consecutive days in different flow cells (FC) on the same sensor chip, with preparation of calibration curve, samples and control samples prior to each analysis (Fig 2). The samples were injected at two dilutions showing a precision of ≤ 2.5% CV within runs and ≤ 4.0% CV between runs. Precision for control samples in buffer runs performed in sextuplicate is also shown in Figure 2.

Control sample	CV (%) within run	CV (%) between runs
High (50 µg/mL)	0.1	0.7
Medium (25 µg/mL)	0.2	1.6
Low (2 µg/mL)	0.4	0.6




Fig 2. Reproducibility of concentration analysis on Sensor Chip Protein A using two sample dilutions and two flow cells including precision for control sample.

## Kinetic characterization example

Tumor necrosis factor alpha (TNF-α) binding kinetics to mAb5 was analyzed. Samples of mAb5 from CHO cell perfusion culture were compared with mAb5 from final fed-batch culture, and purified mAb5, and infliximab (a chimeric monoclonal antibody drug). Perfusion culture samples were collected from 13 d of cell culturing, after cell density had reached steady state. Similar binding kinetics were obtained over the 13 d of perfusion culture. Purification did not significantly alter binding kinetics (Fig 3). All tested antibodies had very high affinity (< 1 nM) with infliximab exhibiting the highest affinity (Fig 3), and binding was very stable with low dissociation rate (Fig 4 and Fig 5). Data was an excellent fit to a 1:1 binding model.




Fig 3. Binding kinetics of TNF-α to mAb5 produced in cell culture, purified mAb5, and infliximab.




Fig 4. Sensorgram (red curve) from analysis of infliximab (ligand) and TNF-α (analyte) and 1:1 fit (black curve).




Fig 5. Sensorgram (red curve) showing analysis of mAb5 from perfusion culture and TNF-α and 1:1 fit (black curve).


## Biacore Sensor Chip Protein G

– a Biacore Extend product



The carboxymethylated dextran matrix on Sensor Chip Protein G is pre-immobilized with a recombinant Protein G—GammaBind™ G, Type 2—which is the same molecule that is used in some protein purification lab resins from GE. The recombinant Protein G binds to a broad range of IgG, such as human (including IgG<sub>3</sub>), rat, rabbit, mouse, guinea pig, goat, sheep, and cow. Antibodies bind only to the heavy chain within the Fc region ensuring antibodies are bound to the surface in a specific orientation.

This is a sensor chip for antibody quantitation and characterization in a wide range of applications.




## Biacore Sensor Chip Protein L

– a Biacore Extend product



The carboxymethylated dextran matrix on Sensor Chip Protein L is pre-immobilized with a recombinant Protein L. Sensor Chip Protein L has the same surface ligand as is used for Capto™ L protein purification products from GE. The recombinant Protein L binds to the human variable domain of kappa light chain subtypes (1, 3, and 4) without interfering with its antigen-binding site. Protein L binds also to kappa light chain from other species such as mouse, rat, and pig (but not to bovine IgG). Sensor Chip Protein L can be used for capture of antibody fragments that lack the Fc region: ScFv, Fab, and Dab.

This is a sensor chip for antibody and antibody fragment quantitation and yes/no screening.



## Biacore Extend products

Biacore Extend products are early access tools providing researchers working in nonregulated environments with the convenience of a wider choice of consumables optimized for Biacore applications. Utilizing the same high-quality components used for the manufacture of standard Biacore consumables, they enable faster and simpler development of a broader range of assays while retaining the high performance levels your work demands.

### What to expect from Biacore Extend products:

- Ordered and supported as standard products
- Instructions for use are found at [gelifsciences.com/instructions](http://gelifsciences.com/instructions)
- QC performed, but no Certificate of Analysis/Conformance issued
- Reduced product information and application examples
- Lot-to-lot variation might be higher compared to standard products
- For some products the delivery time might be up to 8 wk