

Clinical accessories portfolio

## CARESCAPE™ ONE

CARESCAPE ONE Micro-Modules are not clinical accessories.



# ECG

GE Healthcare ECG trunk cables and lead wires help you reliably make the correct connection to GE monitoring devices. These high-quality reusable accessories complement your monitoring equipment, enabling accurate measurements.

ECG cables and leadwires with IEC standard color coding

The list of accessories displayed in this catalogue is not exhaustive, Contact your GE Healthcare representative to obtain the list of all accessories compatible with your equipment.



## ECG LEADWIRES

PRODUCT	CONNECTORS	SIZE	PART NUMBER	QUANTITY
ECG Leadwire set, 3-lead, grabber	Molded	74 cm	2106390-003	1
		130 cm	2106390-004	
ECG Leadwire set, 5-lead, grabber	Molded	74 cm	2106389-003	1
		130 cm	2106389-004	
ECG Leadwire set, 6-lead, grabber	Molded	74 cm	2106398-001	1
		130 cm	2106398-002	
ECG Leadwire set, 4-lead for 10-LD expansion, grabber	Molded	74 cm	2066468-004	1
		130 cm	2066468-008	



2106390-003



2106389-003



2066468-004



2106398-001




**ECG CABLE NEONATALOLOGY**

PRODUCT	CONNECTORS	SIZE	PART NUMBER	QUANTITY
ECG Neonatal Cable 3-LD DIN connector	Molded	3.6 m	2080557-003	1
		1.8 cm	2080557-004	



**ELECTRODES**

PRODUCT	FORM	PART NUMBER	QUANTITY
ECG pre-wired electrode 3-lead, Radiotranslucent - Solid Hydrogel	Square	209101-406	300



**ACCESSORIES**

PRODUCT	PART NUMBER	QUANTITY
DIN Converter, 5-lead for ECG 3/5-lead trunk cables	2106403-001	1
Leadwire separator for ECG Leadwire set, 5-lead	2106399-001	3
Bedsheet clip for ECG trunk cables	2106399-002	1



# TEMPERATURE

Temperature is an important indicator used through the care continuum, from a GP's practice to the ICU. It can indicate infection, response to anesthetic agents, and even tissue perfusion to a limb or specific area. GE Healthcare's temperature probes help you monitor this critical parameter.



## TEMPERATURE PROBES

PRODUCT	PATIENT	PART NUMBER	QUANTITY
Reusable Temperature Probe – Esophageal/Rectal, 280 cm	Adult	2107176-031	1
Reusable Temperature Probe – Esophageal/Rectal, 280 cm	Pediatric	2107176-032	1
Reusable Temperature Probe – Skin, 280 cm	Adult	2107176-033	1

## DISPOSABLE TEMPERATURE PROBES

PRODUCT	PATIENT	PART NUMBER	QUANTITY
Disposable Temperature Probe – General Purpose, 400 series, 9 Fr	Pediatric	M1024229VY	50
Disposable Temperature Probe – General Purpose, 400 series, 12 Fr	Adult	M1024231VY	
Disposable Temperature Probe – Skin, 400 series	Adult - Infant	M1024222VY	



# NON-INVASIVE BLOOD PRESSURE

GE Healthcare provides a comprehensive portfolio of air hoses for use with GE Healthcare monitors connecting to CRITIKON™ blood pressure (BP) cuffs. The air hoses are offered in a variety of lengths and connectors for neonatal and for pediatric/adult use. The hoses are not made with natural rubber latex and help prevent air leakage by providing a secure connection to the blood pressure cuff.

For a complete list of cuffs compatible with your equipment, contact your GE Healthcare sales representative.



## HOSES

AIR HOSES	PATIENT	LENGTH	PART NUMBER	QUANTITY
DINACLICK Connector	Adult - Pediatric	1.2 m	2087389-001	1
		3.6 m	2087389-002	
		7.2 m	2087389-003	
Neo-Snap Connector	Neonate	2.4 m	2089791-002	1
		3,6 m	2089791-003	



2087389-002



2089791-003



## CUFFS

## SOFT-CUF



Our most popular cuff, the CRITIKON™ SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It can be cleaned and reused on multiple patients, but is priced economically so it can be disposed of after a single patient use, helping you affordably prevent infections.

20 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
Infant	Orange	8-13 cm	SFT-P1-2A-5S
Child	Green	12-19 cm	SFT-P2-2A-5S*
Child Long	Green	12-19 cm	SFT-P2-2A-L-5S
Small adult	Royal Blue	17-25 cm	SFT-A1-2A-5S*
Small adult Long	Royal Blue	17-25 cm	SFT-A1-2A-L-5S
Adult	Navy	23-33 cm	SFT-A2-2A-5S*
Adult Long	Navy	23-33 cm	SFT-A2-2A-L-5S*
Large Adult	Pink	31-40 cm	SFT-A3-2A-5S*
Large Adult Long	Pink	31-40 cm	SFT-A3-2A-L-5S
Thigh	Brown	38-50 cm	SFT-T1-2A-5S

Part number with an "-INT" suffix, denoting that the Instructions For Use (IFU) is available in multiple languages.

## ASSORTMENTS DINACLICK CONNECTOR

PRODUCT	2-TUBE DINACLICK
20 Cuffs : 5 Small Adult, 10 Adult, 5 Large Adult	2059303-501
20 Cuffs : 6 Infant, 10 Child, 4 Small Adult	2059303-502
20 Cuffs : 3 Small Adult, 7 Adult, 3 Adult Long, 3 Large Adult, 4 RADIAL-CUF (forearm)	2059303-503
20 Cuffs : 2 Infant, 2 Child, 2 Small Adult, 6 Adult, 2 Adult Long, 2 Large Adult, 2 RADIAL-CUF (forearm), 2 Thigh	2059303-504
4 cuffs: 1 Small Adult, 1 Adult, 1 Large Adult, 1 RADIAL-CUF	2059303-506



CUFFS

# DURA-CUF



CRITIKON™ DURA-CUF™ is gentle and designed for durability. Constructed with materials not made with PVC, making them safer for the environment, and natural rubber latex, making them safer for your patient.

Easy-to-use, fast **DINACLICK connector** enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors. The audible click lets users know there is a proper connection.



5 cuffs/box



SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
Infant	Orange	8-13 cm	DUR-P1-2A-5S
Child	Green	12-19 cm	DUR-P2-2A-5S
Small Adult	Royal Blue	17-25 cm	DUR-A1-2A-5S
Small Adult Long	Royal Blue	17-25 cm	DUR A1-2A-L-5S
Adult	Navy	23-33 cm	DUR-A2-2A-5S
Adult long	Navy	23-33 cm	DUR-A2-2A-L-5S
Large Adult	Wine	31-40 cm	DUR-A3-2A-5S
Large Adult Long	Wine	31-40 cm	DUR-A3-2A-L-5S
Thigh	Brown	38-50 cm	DUR-T1-2A-5S



*The first GE CRITIKON™ blood pressure cuff, DURA-CUF™, combines comfort and durability.*



### ASSORTMENTS DINACLICK CONNECTOR

PRODUCT	2-TUBE DINACLICK
3 cuffs: 1 Small Adult, 1 Adult, 1 Large Adult	2059301-501
6 cuffs : 2 Infant, 3 Child, 1 Small Adult	2059301-502
6 cuffs : 1 Small Adult, 3 Adult, 1 Adult Long, 1 Large Adult	2059301-503
6 cuffs : 1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh	2059301-504



## Validated quality

CRITIKON BP cuffs are validated with GE Healthcare monitors, and are part of an integrated clinical system.



## Better infection control

Source of contamination ignored or underestimated, no other equipment is used more often than a blood pressure cuff without proper cleaning. GE's CRITIKON SOFT-CUF™ cuff provide the optimal solution to protect your patients.



## DINACLICK Connector

Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor faster than other dual-tube connectors.



## Sized for comfort and accuracy

Using the wrong size BP cuff can affect determination accuracy. A full range of CRITIKON™ cuff sizes are available, from neonate, pediatric, adult and obese.

# PRECISE MEASURES FOR INFORMED CHOICES

GE Healthcare offers **the first BP forearm cuff clinically validated<sup>1</sup>** for use on **obese or overly muscular adult patients.**



By adapting to the anatomy of obese patients, the RADIAL-CUF cuff optimizes the accuracy of the **blood pressure** measurements of this particularly risky patient population. **77% of patients** also consider it **more comfortable** than traditional blood pressure cuffs\*.



SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
Obese or overly muscular adult	PURPLE	26-36 CM (FOREARM)	SFT-F1-2A-5S



1. The first clinically validated NIBP Cuff for providing accurate measurements from the forearm, GE Healthcare, 2013.

\* GE Healthcare sponsored and compensated Clinimark LLC for managing and executing the clinical portions of this study (March 2008/ December 2011).





## CUFFS

## NEONATAL

Our neonatal BP cuffs are designed with the performance accuracy you trust and the gentle comfort you expect. Each of our 5 sizes is offered in an extra soft absorbent material (SOFT-CUF™ range) or easy-to-clean vinyl (CLASSIC-CUF™ range). Economical enough for pairing with a single patient throughout hospital stay.

20 cuffs/box			2-TUBE NEO-SNAP	
SIZE	COLOR	LIMB CIRCUMFERENCE	CLASSIC-CUF	SOFT-CUF
1	Orange	3-6 cm	CLA-N1-2B	SFT-N1-2B
2	LT. Blue	4-8 cm	CLA-N2-2B	SFT-N2-2B
3	Green	6-11 cm	CLA-N3-2B	SFT-N3-2B
4	Navy	7-13 cm	CLA-N4-2B	SFT-N4-2B
5	Pink	8-15 cm	CLA-N5-2B	SFT-N5-2B



## ASSORTMENTS NEO-SNAP CONNECTOR


PRODUCT	2-TUBE NEO-SNAP
20 cuffs SOFT-CUF : 2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5	2059305-002
20 cuffs CLASSIC-CUF : 2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5	2059306-002
3 cuffs CLASSIC-CUF : 1 Neo #3, 1 Neo #4, 1 Neo #5	2059306-003



Playful numeric and color-coded animal designs to help you find the right size at a glance.



Actual illustrations may slightly vary.



# SpO<sub>2</sub> SENSORS

Through the legacy of Datex-Ohmeda, GE SpO<sub>2</sub> technology is a story of innovation that has been going on for more than 20 years. TruSignal™ sensors allow continuous non-invasive monitoring of SpO<sub>2</sub> and pulse rate.

We offer a wide range of SpO<sub>2</sub> sensors, designed and tested for their performance and always serving the comfort of your patients.



## TECHNOLOGY

# GE TRUSIGNAL

For your CARESCAPE ONE monitors equipped with GE TruSignal SpO<sub>2</sub> Parameter.


### GE SpO<sub>2</sub> ACCESSORIES

PRODUCT	USE WITH	SIZE	PART NUMBER	QUANTITY
Wide adhesive replacement tape	TS-W-D, TS-SE-3	-	TS-RTW-100	100
Replacement adhesive tape (blue)	TS-AF-10, TS-AF-25	-	TS-RT-100	100
Replacement adhesive tape (bears)			TS-RTB-100	
Foam wrap replacement	TS-SE-3	Small	TS-RWS-24	24
		Medium	TS-RWM-24	
	TS-W-D	Large	TS-RWL-24	
Infant foam sandal (2 to 7,7 kg)	TS-SE-3	-	TS-SND-3	3
Replacement Headband	TS-E-D	-	TS-HB-3	3





REUSABLE SENSORS					
SENSORS	PATIENT	PATIENT WEIGHT	LENGTH	PART NUMBER	QUANTITY
TruSignal Finger sensor	Adult, Pediatric	≥ 20 kg	1 m	TS-F-D	1
TruSignal PediTip sensor	Pediatric	15-30 kg	1 m	TS-SP-D	1
TruSignal FingerTip sensor	Adult	> 30 kg	1 m	TS-SA-D	1
TruSignal Ear sensor	Adult, Pediatric	≥ 10 kg	1 m	TS-E-D	1
TruSignal Wrap sensor	Adult, Pediatric	≥ 3 kg	1 m	TS-W-D	1
TruSignal Sensitive Skin sensor	Adult, Pediatric, Infant, neonate	No limit	1 m	TS-SE-3	3



DISPOSABLE SENSORS					
SENSORS	PATIENT	PATIENT WEIGHT	LENGTH	PART NUMBER	QUANTITY
TruSignal sensor	Adult, Pediatric	≥ 20 kg	0,3 m	TS-AP-10	10
				TS-AP-25	25
TruSignal AllFit sensor	Adult, pediatric, infant, neonate	No limit	0,5 m	TS-AF-10	10
				TS-AF-25	25
TruSignal Adhesive Wrap sensor	Adult	> 20 kg	0,3 m	TS-AAW-10	10
				TS-AAW-25	25
TruSignal Adhesive Wrap sensor	Pediatric, infant	3 - 20 kg	0,3 m	TS-PAW-10	10
				TS-PAW-25	25



## TECHNOLOGY

## NELLCOR™


For your CARESCAPE One monitors equipped with with Nellcor SpO<sub>2</sub> Parameter.



## REUSABLE NELLCOR SENSORS

SENSORS	MANUFACTURER REFERENCE	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
Nellcor DuraSensor Finger Probe	DS 100A	Adult	< 30 kg	70124021	1
Nellcor SpO <sub>2</sub> Multisite Sensor	D-YS	Adult	> 1 kg	70124033	1
Nellcor SpO <sub>2</sub> Sensor	OXI-A/N	Neonate	< 3 kg	70124035	1
		Adult	> 40 kg		
Nellcor SpO <sub>2</sub> Sensor	OXI-P/I	Pediatric, infant	3-40 kg	414248-001	1 delivered with 50 adhesive strips

## DISPOSABLE NELLCOR SENSORS



SENSORS	MANUFACTURER REFERENCE	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
Nellcor SpO <sub>2</sub> Finger Sensor	MAX-A	Adult	> 30 kg	70124027	24
Nellcor SpO <sub>2</sub> Sensor	MAX-N	Adult	> 40 kg	70124032	24
		Neonate	< 3 kg		

## FIXATION CLIP



PRODUCT	PART NUMBER	QUANTITY
Nellcor SpO <sub>2</sub> Ear-Clip D-YSE, for use with 70124033	70124034	1

## ADHESIVE TAPE



PRODUCT	PART NUMBER	QUANTITY
Nellcor adhesive tape (ADH-A/N) - Ref: 70124035 (OXI-A/N)	2016130-001	100

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.




## TECHNOLOGY

# MASIMO™


For your CARESCAPE One monitors equipped with MASIMO SpO<sub>2</sub> Parameter.



## DISPOSABLE SENSORS



2027253-001



SENSORS	MANUFACTURER REFERENCE	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
Adhesive LNCS sensors	LNCS NEO	Neonate	< 3 kg	2051069-001	20
	LNCS NEOPT	Neonate	< 1 kg	2051070-001	
	LNCS Adtx	Adult	> 30 kg	2027253-001	
Digit RD SET sensor	RD SET Adt	Adult	> 30 kg	MASIMO-RD 4000	20
	RD Set Inf	Infant	3-20 Kg	MASIMO-RD 4002	
	RD Set Neo	Neonate/Adult	< 3 or > 40 kg	MASIMO-RD 4003	
	RD SET NeoPt	Neonate	< 1 kg	MASIMO-RD 4004	
Digit RD Rainbow SET Sensor	RD Rainbow Set -2 Adt	Adult	> 30 kg	MASIMO-RD 4026	20
	RD Rainbow Set -2 Pdt	Pediatric	10-50 kg	MASIMO-RD 4027	
	RD Rainbow Set -2 Inf	Infant	3-20 Kg	MASIMO-RD 4028	
	RD Rainbow Set -2 Neo	Neonate/Adult	< 3 or > 40 kg	MASIMO-RD 4029	

## REUSABLE SENSORS



Masimo-RD-4051

SENSORS	MANUFACTURER REFERENCE	PATIENT	PART NUMBER	QUANTITY
MASIMO RD SET	RD SET DCI	Adult	MASIMO-RD-4050	1
	RD SET DCI-P	Pediatric	MASIMO-RD-4051	

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.



**REUSABLE LNCS SENSORS**

SENSORS	MANUFACTURER REFERENCE	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
Masimo Finger sensor	LNCS DC-I	Adult	> 30 kg	2027258-001	1
	LNCS DC-IP	Pediatric	10-50 kg	2027259-001	
Masimo TipClip Ear sensor	LNCS TC-I	Adult	> 30 kg	2027261-001	1



**REUSABLE LNCS CABLES**

CABLES	MANUFACTURER REFERENCE	LENGTH	PART NUMBER	QUANTITY
MASIMO LNCS CABLE	RED 25 LNC-10	3 m	2106257-001	1
	RED 25 LNC-4	1.2 m	2106257-002	
	RED 25 LNC-1	30 cm	2106257-003	




**MASIMO RD CABLE**

CABLE	MANUFACTURER REFERENCE	LENGTH	PART NUMBER	QUANTITY
MASIMO RD CABLE	RD Rainbow SET 25R-1.5	0.45 m	MASIMO-RD-4074	1
	RD Rainbow SET 25R-05	1.5 m	MASIMO-RD-4076	
	RD Rainbow SET 25R-08	2.4 m	MASIMO-RD-4106	
	RD Rainbow SET 25R-12	3.6 m	MASIMO-RD-4078	



These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Medtronic  
Microstream CO<sub>2</sub>



# CARESCAPE CO<sub>2</sub>- MICROSTREAM™

CARESCAPE CO<sub>2</sub>-Microstream parameter for measuring end-tidal carbon dioxide (EtCO<sub>2</sub>), FiCO<sub>2</sub>, and respiration rate with Medtronic Microstream™ technology.

## INTUBATED PATIENT

DESCRIPTION	PATIENT	LENGTH	PART NUMBER	QUANTITY
Microstream™ Advance Intubated CO <sub>2</sub> Filter Line, Extended Duration	Adult-Pediatric	200 cm	MVAIH	25
Microstream™ Advance Intubated CO <sub>2</sub> Filter Line, High Humidity	Adult-Pediatric		MVAIHH	
Microstream™ Advance Intubated CO <sub>2</sub> Filter Line, Extended Duration	Neonatal-Infant		MVIIH	
Microstream™ Advance Intubated CO <sub>2</sub> Filter Line, High Humidity	Neonatal-Infant		MVIIHH	

MVAIH, MVAIHH: Extended drier moisture reduction technology for high ambient humidity.

MVIIH, MVIIHH: Intended for use with ETT > 4.5mm. Extended drier moisture reduction technology for high ambient humidity.

These accessories are sold ONLY with the equipment.



**INTUBATED PATIENT, SHORT-TERM USE**

DESCRIPTION	PATIENT	LENGTH	PART NUMBER	QUANTITY
Microstream™ Advance Intubated CO <sub>2</sub> Filter Line, Short-term use: Procedural/Emergency	Adult-Pediatric	200 cm	MVAI	25

Intended for use with ETT size > 4.5mm

**NON INTUBATED PATIENT, SHORT-TERM USE**

DESCRIPTION	PATIENT	LENGTH O <sub>2</sub> TUBING FEMALE CONNECTOR	PART NUMBER	QUANTITY
Microstream™ Advance Oral-Nasal CO <sub>2</sub> Filter Line with O <sub>2</sub> Tubing	Adult	200 cm	MVAO	25
	Adult	400 cm	MVAOL	
Microstream™ Advance Oral-Nasal CO <sub>2</sub> Filter Line with O <sub>2</sub> Tubing	Pediatric	200 cm	MVPO	


Samples from both the mouth and the nose. Integrated O<sub>2</sub> tubing for supplemental oxygen delivery

**NON INTUBATED PATIENT, EXTENDED DURATION**

DESCRIPTION	PATIENT	LENGTH O <sub>2</sub> TUBING FEMALE CONNECTOR	PART NUMBER	QUANTITY
Microstream™ Advance Oral-Nasal CO <sub>2</sub> Filter Line with O <sub>2</sub> Tubing	Adult	200	MVAOH	25
	Neonatal	300	MVINOH	

Samples from both the mouth and the nose. Drier moisture reduction technology. Integrated O<sub>2</sub> tubing for supplemental. Oxygen delivery





# CARESCAPE rSO<sub>2</sub>-INVOS™

CARESCAPE rSO<sub>2</sub> - INVOS parameter for measuring regional oxygen saturation of blood (rSO<sub>2</sub>) in cerebral and somatic tissues.



## INVOS SENSOR

DESCRIPTION	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
INVOS™ Cerebral/Somatic Oximetry Sensors	Adult	> 40 kg	PMSENS71-A-10	10

## INVOS CABLE

DESCRIPTION	PATIENT	PART NUMBER	QUANTITY
INVOS™ Reusable Sensor Cable	Adult	PMAC71RSC	1

These accessories are sold ONLY with the equipment.



# CO<sub>2</sub> ACCESSORIES

We offer a range of Airway Adapter sets dedicated in LoFlo and help with humidity problem, they are designed and tested for their performance and always serving the comfort of your patients.



## AIRWAY ADAPTER

DESCRIPTION	PATIENT TYPE	LENGTH	PART NUMBER	QUANTITY
Kit LoFlo Airway Adapter & Sample Line	Adult/Pediatric	2 m	2013068-001	10
	Pediatric/Infant		2013068-002	
Kit LoFlo Airway Adapter & Sample Line with Dehumidification Tubing	Adult/Pediatric		2013068-003	
	Pediatric/Infant		2013068-004	

2013068-001



2013068-002



2013068-004





# INVASIVE BLOOD PRESSURE

GE Healthcare offers compatibility with four different invasive blood pressure transducers by partnering with Argon, Edwards Lifesciences TruWave, ICU Medical Transpac-IV and Utah/Biotrans. When your CARESCAPE™ ONE is equipped with the Generic Interface micro-module (2097273-009)\* you can continue to utilize the legacy GE Healthcare invasive blood pressure cables you've depended on for years. We offer generic interface cables for these four transducer manufacturers in two different lengths to fit your needs.



## INVASIVE BLOOD PRESSURE CABLES



PRODUCT	TYPE	LENGTH	PART NUMBER	QUANTITY
Invasive Blood Pressure cable - Utah/Biotrans	Single	3,6 m	2104107-001	1
		1.2 m	2104107-002	
Invasive Blood Pressure cable - ICU Medical Transpac-IV	Single	3,6 m	2104158-001	1
		1.2 m	2104158-002	
Invasive Blood Pressure cable - Edwards LS TruWave	Single	3,6 m	2104162-001	1
		1.2 m	2104162-002	
Invasive Blood Pressure cable - Argon	Single	3,6 m	2104166-001	1
		1.2 m	2104166-002	

\*The above cables will require the Generic interface cable 2097273-009 for CARESCAPE ONE.



## About GE Healthcare

GE Healthcare is a leading global medical technology and digital solutions innovator, GE Healthcare enables clinicians to make faster, more informed decisions through intelligent devices, data analytics, applications and services, supported by its Edison intelligence platform. With over 100 years of healthcare industry experience and around 50,000 employees globally, the company operates at the center of an ecosystem working toward precision health, digitizing healthcare, helping drive productivity and improve outcomes for patients, providers, health systems and researchers around the world. Follow us on Facebook, LinkedIn, Twitter and Insights , or visit our website [www.gehealthcare.com](http://www.gehealthcare.com) for more information.

The list of accessories displayed in this catalogue is not exhaustive, Contact your GE Healthcare representative to obtain the list of all accessories compatible with your equipment.

[www.gehealthcare.co.uk](http://www.gehealthcare.co.uk)

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information.

Data subject to change.

GE, the GE Monogram, CARESCAPE, CRITIKON, DURA-CUF, SOFT-CUF, DINACLICK and TruSignal are trademarks of General Electric Company. Microstream and Invos are trademarks of Medtronic.

\*All other third-party trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE.

Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

© 2021 General Electric Company. All rights reserved.

JB00228XE Rev4

