



Giraffe



Clinical accessories portfolio

## Infant Care

[gehealthcare.co.uk](http://gehealthcare.co.uk)


Not available in all countries



Maintaining a clean airflow and a dark, quiet environment is important for neonatal patients. Equally important is keeping the baby's temperature stable.

ACCESSORIES		
DESCRIPTION	PART NUMBER	QUANTITY
Air Filters*	2107356-001	10
Bag and Mask Holder Clip	2107357-001	1

\*Compatible with incubators only.



### GIRAFFE™ HOOD COVERS

DESCRIPTION	PART NUMBER	QUANTITY
OmniBed™ Cover	6600-0825-800	1
Incubator Cover	6600-0846-800	1

Sold with equipment only in some countries.



### TEMPERATURE PROBES AND PATCHES

DESCRIPTION	PART NUMBER	QUANTITY
Giraffe and Panda Disposable Probes	2074816-001	10
Giraffe and Panda Disposable Probes	2074817-001	50
Heat Reflecting Patches	0203-1980-300	50

Compatible with Giraffe Warmer, Giraffe Omnibed, Giraffe Incubator Carestation, and Panda Warmer. Sold with equipment only in some countries.



0203-1980-300

Not available in all countries


The BiliSoft Pad Covers and Pad Nests are designed for use with the BiliSoft Phototherapy Systems, which are used to treat both premature and full-term infants for hyperbilirubinemia.

### BILISOFT™ DISPOSABLE PAD COVERS

DESCRIPTION	BILISOFT™	BILISOFT™ 2.0	SIZE	PART NUMBER	QUANTITY
Disposable Pad Covers, Flat		X	Small	2105422-001	20
			Small	2105418-001	50
	X		Small	2107201-001	20
	X		Small	2107200-001	50
	X	X	Large	2105423-001	20
	X	X	Large	2105419-001	50
Disposable Nest Pad Covers		X	Small	2105420-001	15
	X		Small	2107202-001	
	X	X	Large	2105421-001	



Flat cushioned cover with straps  
The flat, cushioned BiliSoft™ Cover lets you swaddle baby and phototherapy pad together.



Cushioned cover with developmental positioning foot roll and straps  
The BiliSoft™ Nest offers sick babies the boundaries and support they need.

### Product compatibility

The BiliSoft™ Pad Covers and Nest are designed for use with the BiliSoft™ and BiliSoft™ 2.0 Phototherapy System, which is used to treat both premature and full-term infants for hyperbilirubinemia.

Available in small and large sizes to match the size of the Light Pad\*

\* Note: the BiliSoft™ 2.0 Disposable Pad Cover (Flat and Nest) in Size Large also fits the BiliSoft™ Phototherapy System (Legacy) large pad.

Specifications available on request.



BiliSoft™

BiliSoft™ 2.0

Not available in all countries



**T-Piece Circuit: Disposable resuscitator designed for neonates, provides accurate reliable ventilation and pulmonary resuscitation.**



2109061-001  
Disposable T Piece with Mask

**INFANT RESUSCITATOR T-PIECE CIRCUIT KIT**

DESCRIPTION	MASK SIZE	PART NUMBER	QUANTITY
Disposable T-Piece	No mask	2109059-001	20
Disposable T-Piece with Mask	0	2109060-001	10
	1	2109061-001	10

**Product compatibility**

Compatible with GE Healthcare resuscitation devices using a standard connector.

Specifications available on request.



Panda™ Warmer



Giraffe™ Stand-alone Infant Resuscitation System

Not available in all countries

**Providing a stable, thermoregulated environment, the Lullaby™ Warmer helps clinicians meet the requirements of different care areas.**

**TEMPERATURE PROBES**

DESCRIPTION	PART NUMBER	QUANTITY
Disposable Temperature Probe	6600-0196-700	50
	6600-0208-700	10
Reusable Temperature Probe – Lullaby™ Warmer	2044384-001	1

Sold with equipment only in most countries.



6600-0196-700

**TEMPERATURE PROBE PATCHES**

DESCRIPTION	PART NUMBER	QUANTITY
Heat Reflecting Patches	0203-1980-300	50

Sold with equipment only in most countries.



0203-1980-300

Not available in all countries





Providing big benefits for your tiniest patients, GE Healthcare's neonatal cuffs are designed with the accuracy you trust and the gentle comfort you expect. We offer a full range of neonatal sizes with the choice of extra soft, absorbent material or easy-to-clean vinyl. The playful, color-coded animal designs help you find the right size at a glance and they are economically priced for single-patient-use.



Neonatal			2-Tube Neo-Snap Connector		1-Tube Neo-Snap Connector*	
Size	Color	Limb Circumference	SOFT-CUF	CLASSIC-CUF	SOFT-CUF	CLASSIC-CUF
1	Orange	3-6 cm	SFT-N1-2B	CLA-N1-2B	SFT-N1-1B	CLA-N1-1B
2	Light Blue	4-8 cm	SFT-N2-2B	CLA-N2-2B	SFT-N2-1B	CLA-N2-1B
3	Green	6-11 cm	SFT-N3-2B	CLA-N3-2B	SFT-N3-1B	CLA-N3-1B
4	Navy	7-13 cm	SFT-N4-2B	CLA-N4-2B	SFT-N4-1B	CLA-N4-1B
5	Rose	8-15 cm	SFT-N5-2B	CLA-N5-2B	SFT-N5-1B	CLA-N5-1B

NEONATAL ASSORTMENT PACK	CUFF TYPE	PART NUMBER	DESCRIPTION
2-Tube Neo-Snap Connectors	SOFT-CUF	2059305-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-003	1 Neo #3, 1 Neo #4, 1 Neo #5 (3 cuffs/pack)
1-Tube Neo-Snap Connectors*	SOFT-CUF	2059305-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)
	CLASSIC-CUF	2059306-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/box)

Note: 20 cuffs/box

\* 1-Tube Neo-Snap connectors cannot be used on GE monitoring equipment.

Not available in all countries

GE Healthcare has a wide variety of monitoring solutions, with a full compliment of products, for all neonates and infants. For your CARESCAPE B450 / B650 / B850 monitors equipped with a PDM module with MASIMO SpO<sub>2</sub> technology.



MASIMO™ LNCs & RD SET INTERCONNECT CABLES

DESCRIPTION	SIZE	MASIMO™ REFERENCE	PART NUMBER	QUANTITY	
LNC-10 Interconnect Cable	3 m/10 ft	For GE DINAMAP	1814	2027263-001	1
		For GE	2016	2027263-002	
RD SET Interconnect Cable	3.6 m/12 ft	4085	MASIMO-RD-4085	1	

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.



2051070-001

MASIMO™ LNCs & RD SET DISPOSABLE SENSORS

DESCRIPTION	PATIENT WEIGHT	MASIMO™ REFERENCE	PART NUMBER	QUANTITY	
LNCs SpO <sub>2</sub>	Infant	3-20 kg	2328	2051068-001	20
	Neonatal PT	< 1 kg	2330	2051070-001	20
	Adult/ Neonatal	> 40 kg/ < 3 kg	2329	2051069-001	20
RD SET® SpO <sub>2</sub>	Infant	3-20 kg	4002	MASIMO-RD-4002	20
	Adult/ Neonatal	> 40 kg/ < 3 kg	4003	MASIMO-RD-4003	20
	Premature Baby	< 1 kg	4004	MASIMO-RD-4004	20

Requires Compatible Interconnect Cable.

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

MASIMO™ LNCs & RD SET SENSOR ACCESSORIES

DESCRIPTION	MASIMO™ REFERENCE	PART NUMBER	QUANTITY	
LNCs Replacement Tape	Infant	2307	2051071-001	100
	Neonatal	2308	2051073-001	100
	Neonatal PT	2309	2051074-001	100
RD Set Replacement Tape	Infant	4046	MASIMO-RD-4046	102
	Neonatal	4047	MASIMO-RD-4047	102
RD Set Replacement Wrap	4048	MASIMO-RD-4048	10	


These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices in some countries. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Not available in all countries



For CARESCAPE™ B450/B650/B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO<sub>2</sub> technology.

NELLCOR™ INTERCONNECT CABLES				
Description	CONNECTOR	LENGTH	PART NUMBER	QUANTITY
OxiMax™ SpO <sub>2</sub> Interconnect Cable	GE	1.2 m/4 ft	2021406-002	1
		3 m/10 ft	2021406-001	1
Interconnect Cables*	DOC-4	1.2 m/4 ft	2059000-001	1
	DOC-8	2.4 m/8 ft	2059001-001	1
	DOC-10	3 m/10ft	2008773-001	1



\*These accessories are sold ONLY with the equipment in some countries.

NELLCOR™ REUSABLE SENSORS					
SENSORS	PATIENT	PATIENT WEIGHT	NELLCOR™ REFERENCE	PART NUMBER	QUANTITY
Two-piece SpO <sub>2</sub> sensor	Neonatal/Adult	< 3 kg/> 40 kg	OXI-A/N	70124035	1 (delivered with 24 adhesive strips)
	Pediatric, infant	3-40 kg	OXI-P/I	414248-001	1 (delivered with 24 adhesive strips)


Requires Compatible Interconnect Cable.  
These accessories are sold ONLY with the equipment.



Not available in all countries

For CARESCAPE™ B450/B650/B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO<sub>2</sub> technology.

NELLCOR™ DISPOSABLE SENSORS					
SENSORS	NELLCOR™ REFERENCE	PATIENT	PATIENT WEIGHT	PART NUMBER	QUANTITY
Adhesive SpO <sub>2</sub> sensor	MAXP	Pediatric	10-50 kg	70124022	24
	MAXI	Infant	3-20 kg	70124026	24
	MAX-N	Adult	> 40 kg	70124032	24
Neonate		< 3 kg			



Requires Compatible Interconnect Cable.  
These accessories are sold ONLY with the equipment.

NELLCOR™ SENSOR ADHESIVE TAPE					
PRODUCT	PATIENT	NELLCOR™ REFERENCE	USE WITH	PART NUMBER	QUANTITY
Sensor Wrap Replacement Tape	Adult/Neonatal	ADH-A/N	OXI-A/N, D-YS	2016130-001	100
	Pediatric/Infant	ADH-P/I	OXI-P/I, D-YS	2016131-001	100




These accessories are sold ONLY with the equipment.

Not available in all countries



ECG CABLE			
DESCRIPTION	LENGTH	PART NUMBER	QUANTITY
3-lead ECG trunk cable Neonatal DIN	1,2 m	2106306-004	1
	3,6 m	2106306-003	



ELECTRODES			
DESCRIPTION	FORM	PART NUMBER	QUANTITY
ECG pre-wired electrode 3-lead, Radiotranslucent - Solid Hydrogel	Square	2009101-406	300

**About GE Healthcare**

GE Healthcare is a leading global medical technology and digital solutions innovator, GE Healthcare enables clinicians to make faster, more informed decisions through intelligent devices, data analytics, applications and services, supported by its Edison intelligence platform. With over 100 years of healthcare industry experience and around 50,000 employees globally, the company operates at the center of an ecosystem working toward precision health, digitizing healthcare, helping drive productivity and improve outcomes for patients, providers, health systems and researchers around the world. Follow us on Facebook, LinkedIn, Twitter and Insights, or visit our website [www.gehealthcare.com](http://www.gehealthcare.com) for more information.

The list of accessories displayed in this catalogue is not exhaustive. Contact your GE Healthcare representative to obtain the list of all accessories compatible with your equipment.

[www.gehealthcare.co.uk](http://www.gehealthcare.co.uk).

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information.

Data subject to change.

GE, the GE Monogram, Giraffe, Panda, Lullaby, and Giraffe Shuttle are trademarks of General Electric Company.

\*All other third-party trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE.

Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

© 2022 GE. All rights reserved.

JB04900XE